## PRODUCT CODES

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## ENGLISH

**Description**  
 Reusable label that reminds the user that a period of time has elapsed.

- Instructions**
- A** Set period Timertag will measure.
  - B** Green light.
  - C** Button to start, stop and test.
  - D** Red light.

Fix magnetically, or using removable adhesive.

**Start** Press and hold button until green light flashes three times.

**Stop** Press and hold button until red light flashes three times.

- Test** Press button <1 second:
- Green blink (<1s) = Ready.
  - Green light (>1s) = Active.
  - Red light = Expired.
  - No lights = Expired.

**Alert (Timertags 2-72h)**  
 Red light flashes every second for 2h, then once every 10s for 22h.

**Alert (Timertags 7d-1y)**  
 Red light flashes every second for 8h, then once every 10s for 14d.

**Cautions**  
 Only use dry-erase (whiteboard) pens on Timertags.

**Disposal**  
 Dispose via normal e-waste recycling. RoHS compliant in accordance with EU Directive 2011/65/EU.

**Warnings**  
 Avoid X-Rays or Magnetic Resonance Imaging (MRI). Contains a lithium manganese battery. Seek medical attention if swallowed. Keep out of reach of children. Do not burn, heat above 100°C (212°F), puncture or crush - explosion risk.

### Symbols

- Location.
- Date and time.
- Insertion.
- Reason.
- Manufacturer.
- Expiry date.
- Batch code.
- Authorised representative in EU.
- Do not burn or incinerate.
- Contains lithium battery.
- Designed to operate between 10°C (50°F) and 50°C (122°F).
- Keep out of reach of children.
- Instructions for use.
- Do not heat above 100°C (212°F); explosion risk.
- Do not dispose in normal waste.
- Do not use if package damaged.

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