
## PRODUCT CODES

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## EN

**Description** Single-use label alert that reminds the user that a period of time has elapsed.

**Illustration - A** Set period Timertag will measure. **B** Activation tag. **Start** Pull out activation tag. Green light will flash once every second for 30 seconds, then once every 5 minutes until the set time period has been reached.

**Place** Remove the release liner from the rear of the Timertag and affix to a smooth laminated or polished surface.

**Alert** After reaching the set time period, a red light will flash twice every three seconds for 24 hours, then once every 5 seconds for at least 7 days.

**Cautions** Water resistant only. Not waterproof. Do not immerse. Adhesive is removable from polished or laminated surfaces only. Do not place on painted or textured surfaces.

**Disposal** Remove Timertag by pulling on the recessed edge and twisting away from the underlying surface. Dispose via normal e-waste recycling. RoHS compliant in accordance with EU Directive 2011/65/EU.

**Cleaning** Timertags are compatible with surface wipes containing: Detergent, Hypochlorite (1%), Peracetic acid, Isopropyl alcohol, Ethyl alcohol, Chlorhexidine, Quaternary Ammonium Compounds, Hydrogen peroxide, 2- Butoxyethanol, Polyhexamethylene biguanide (PHMB), Triamine, Isopropanolol, Citric Acid, Chlorine dioxide, Iodine.

**Shelf Life** 5 years from date of manufacture.

**Performance** Median error (ME)=0.5%. Maximum observed error (MOE)=1.5%. TTBM90D ME=0.45 days MOE=1.35 days. TTBM180D ME=0.9 days MOE=2.7 days. TTBM365D ME=1.8 days MOE=5.4 days.

**Warnings** Avoid X-Rays or Magnetic Resonance Imaging (MRI). Contains a lithium manganese battery. Seek medical attention if swallowed. Keep out of reach of children. Do not burn, heat above 100°C (212°F), puncture or crush - explosion risk. Timertags are workflow aides and do not prevent disease - no health benefit claims are made.

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