

# The Guardian CPV

## Laryngeal Mask - The Ultimate Seal

### Features & Benefits



- **Vented mask** – airway access (std 15mm connection) and gastric access (up to 16Ga NG tube)
- **All silicone cuff and shaft** – softer / flexible / warms quickly insitu / better fit / allows head movement without losing seal
- **Pre-shaped** – for easier insertion
- **Airway lumen fits an Aintree** – able to secure airway in an emergency
- **Three dimensional cuff inflation** – inflates anterior and posterior for better seal of up to 40cmH2O
- **Integrated bite block** – maintains patent airway
- **Cuff pilot valve** – permanent manometer for accurate cuff pressures
- **Scavenge port** – quick wall suction access for suctioning at tip of cuff
- **No tab** – less risk of pressure area peri lip area
- **Single use** – no delays in recycling
- **Price competitive** – savings of approximately 20% on existing disposable vented mask options

Distributor: Queensland / Northern NSW



MEDICAL TECHNOLOGY AND SERVICES

Phone: + 617 3139 1207

Fax: + 617 3036 6187

Mobile: 0401 422 767

Email: [sales@lateralmedical.com](mailto:sales@lateralmedical.com)