


## Benninger Microlaryngeal Instruments

### A New Level of Versatility

Laryngeal instruments that combine quality, precision, and stability

The Benninger Microlaryngeal Instruments are designed to provide maximum stabilisation of the surgeon's hand. The reverse action handles are thumb-activated to reduce the amount of distal tip movement and wrist movement required by the surgeon for easier, atraumatic dissection.



This set includes forceps, knives, retractors, injection instruments, suction cannulas, and spatulas; everything required for microlaryngeal manipulation and dissection

### ordering information

**MCLB2**      **MicroFrance® Benninger Microlaryngeal Instrument Set**      1 set

*Set Includes:*

MCL-S21R	Grasper, Right, Heart-Shaped	1
MCL-S22R	Grasper, Left, Heart-Shaped	1
MCL-S23R	Alligator Forceps, Right	1
MCL-S24R	Alligator Forceps, Left	1
MCLS18R	Cup Forceps, Downbiting, 2.3 mm	1
MCLB94	Laryngeal Needle Driver, Angled 70° Up, Ratchet Lock	1
MCLB30	Laryngeal Suture Knot Pusher	1
MCLB16	Laryngeal Micro-Scissors, Curved Right, 1.5 mm Blades	1
MCLB17	Laryngeal Micro-Scissors, Curved Left, 1.5 mm Blades	1
MCLB31	Laryngeal Micro-Scissors, Angled 70° Up, 1.5 mm Blades	1

**Sataloff Microlaryngeal Instruments Included in this Set:**

MCL-S0	Long Connector for Injection Needle	1
MCL-S1	Spatula, Curved, 1 mm x 9 mm	1
MCL-S7	Spatula, Straight, 1 mm x 10 mm	1
MCL-S4	Injection Needle, Straight	6
MCL-S25	Sickle Knife (single-use) For use with MCL-S19	6
MCL-S19	Universal Instrument Handle, Ring Handle with Finger Loop	1
MCL-S16	Suction, 3 French	1
3734109	Suction, 6 French	1
MCLS30	Knife Handle	2
MCLB2SC	Instrument Tray	1



**ENT**

MICROFRANCE INSTRUMENTS