

## Bouchayer Laryngoscopes

A properly designed laryngoscope is a necessity for effective direct laryngoscopy. The Bouchayer Laryngoscopes offer a large enough proximal opening for true binocular microscope vision, an elevated distal tip for retraction of the epiglottis and proper exposure of the anterior commissure, and an integrated light fiber for optimal visualisation.

### ordering information

**MCL700 Bouchayer Laryngoscope Set, Large** 1

Set includes:

MCL77 Bouchayer Laryngoscope, Large 1

MCL78 Light Fiber 1

**MCL710 Bouchayer Laryngoscope Set, Medium** 1

Set includes:

MCL79 Bouchayer Laryngoscope, Medium 1

MCL78 Light Fiber 1

### Bouchayer Laryngoscope Accessories

MCL76 Laryngoscope Suspension Apparatus 1

Overall Length 35 cm (13.8")

MCL35 Chest Plate for Suspension Apparatus 1


60 mm Wide x 140 mm Long

MCL78 Replacement Light Fiber for MCL77 and MCL79 1

3717023 Laryngoscope Instrument Tray 1

52 .5 cm L x 25 cm W x 17 cm D

*Holds Two Laryngoscopes, Two Fiberoptic Cables, Suspension Apparatus, and Chest Plate*



**ENT**

MICROFRANCE INSTRUMENTS