

Sataloff Laryngoscopes

MicroFrance® Sataloff Laryngoscopes enable optimal visualization with an enhanced design for examining the vocal cords at the anterior commissure and vallecular regions.

A large proximal opening allows true binocular vision through a microscope, while dual light fibers and dual suction offer maximum illumination and visualization during microsurgery. The elevated triangular tips are blunted to reduce tissue trauma. All laryngoscopes can be used with the MicroFrance suspension apparatus and chest plate.

ordering information

MCL120	Laryngoscope, Pediatric	1
MCL121	Laryngoscope, Medium Male	1
MCL123	Laryngoscope, Medium Female	1
MCL122	Laryngoscope, Large Adult	1
MCL124	Vallecula Laryngoscope	1
MCL125	Anterior Commissure Laryngoscope	1

For patients with limited anterior exposure

MCL76	Suspension Apparatus	1
MCL35	Chest Plate for Suspension Apparatus	1
1853100	Fiberoptic Cable	1
3717023	Laryngoscope Instrument Tray 52.5 cm L x 25.0 cm W x 17.0 cm D	1

Holds two laryngoscopes, two fiberoptic cables, suspension apparatus, and chest plate

Anterior Commissure Laryngoscopes

MCL120, MCL121, MCL122, MCL123

- For the majority of patients
- Available in a range of sizes

MCL125

- For patients with limited laryngeal exposure
- Shorter height on the proximal opening and more acute angle on the distal tip

Vallecula Laryngoscope

MCL124

- Wide proximal opening and high blunted lift on the distal tip
- Provides ample working space and visibility within the vallecular region



ENT

MICROFRANCE INSTRUMENTS