

## Terris Thyroidectomy Instrument Set

The Terris Thyroidectomy Instrument Set offers an innovative approach to address the unique challenges of thyroid surgery. These instruments may be used for both conventional “open” and minimally invasive techniques.

Angled, flared retractors provide maximum access and visibility without interfering with the anterior chest wall

Malleable suction elevator enables tissue retraction and ultrasonic steam evacuation

Malleable blunt tip suction is customised for atraumatic suction around delicate nerves

Low profile elevators are designed for precise dissection of recurrent laryngeal nerve

Slender peanut holder design allows for unobtrusive blunt dissection; one-handed mechanism enables rapid loading and unloading



### ordering information

MCLT00	Thyroidectomy Instrument Set	1 set
MCLT10	Retractor, Small	1 ea
MCLT11	Retractor, Large	1 ea
MCLT20	Elevator, 2mm	1 ea
MCLT21	Elevator, 4mm	1 ea
MCLT30	Peanut Holder	1 ea
MCLT40	Malleable Suction Elevator	1 ea
MCLT50	Malleable Nerve Suction	1 ea
3717012	Medium Instrument Tray	1 ea



**ENT**

MICROFRANCE INSTRUMENTS