

Vaughan Suction Sinus Instruments

Visibility is critically important in sinus surgery, and blood in the operative field is a significant problem. The innovative MicroFrance® Vaughan Suction Sinus Instruments evacuate blood from the surgical site for optimal visibility in the most challenging areas of the sinuses, such as the sphenoid sinus and frontal recess.

The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

- Suction mushroom punches are designed with the suction at the tip of the instrument
- A unique cleaning tool allows soft tissue to be removed quickly and safely during surgery
- Suction picks aid in the dissection of closed cells while evacuating blood
- Malleable suctions can be bent and shaped to fit challenging anatomical constraints

ordering information

MCEN116	Suction Mushroom Punch, 4.0 mm, Straight	1
MCEN118	Suction Mushroom Punch, 6.0 mm, Straight	1
MCEN117	Suction Mushroom Punch, 4.0 mm, 60o Angled Up	1
MCEN110	Mushroom Punch Cleaner	1
MCEN120	Suction Sickle Knife	1
MCEN111	Suction Pick, 90o	1
MCEN112	Suction Pick, 45o Angled Back	1
MCEN113	Suction Ring Curette	1
MCEN114	Suction Cup Forceps	1
MCEN115-6	Malleable Suction, 6 French	1
MCEN115-8	Malleable Suction, 8 French	1
MCEN115-9	Malleable Suction, 9 French	1



ENT

MICROFRANCE INSTRUMENTS