

Wormald Bipolar Forceps with Suction

Available in two lengths, the MicroFrance® Wormald Bipolar Forceps with Suction help achieve hemostasis during endoscopic sinus and anterior skull base surgery. The longer forceps are able to reach the pituitary and come in straight and up-curved versions. These forceps feature a slim profile and fine tips to provide optimum access within difficult-to-reach areas. The integrated suction can be advanced independently of the cautery tips to provide clear visibility of the operative site, and then retracted from the line of sight during cauterization. An integrated flush port facilitates cleaning and debris removal.



ordering information

MCEN64K. Wormald Bipolar Forceps with Suction Set, 105 mm. 1 Set

Set includes:

MCEN64	Wormald Bipolar Forceps with Suction, 105 mm	1
MCEN640	Instrument Tray	1
CP393	Bifurcated (Banana Plug) Bipolar Cable	1

MCEN641K. Wormald Bipolar Forceps with Suction Set, 170 mm. 1 Set

Set includes:

MCEN641	Wormald Bipolar Forceps with Suction, 170 mm	1
MCEN640	Instrument Tray	1
CP393	Bifurcated (Banana Plug) Bipolar Cable	1

MCEN643K. Wormald Bipolar Forceps with Suction Set Up-curved, 170 mm. 1 Set

Set includes:

MCEN643	Wormald Bipolar Forceps, Up-curved, 170 mm	1
MCEN644	Disposable Suction Tip	5
MCEN640	Instrument Tray	1
CP393	Bifurcated (Banana Plug) Bipolar Cable	1



ENT

MICROFRANCE INSTRUMENTS