

LT-300 HD

High Definition Digital Video Colposcope

Features

- ⌘ 3.27 mega pixels
- ⌘ SONY digital camera
- ⌘ Mini LCD viewing screen
- ⌘ Adjustable LED lights
- ⌘ Progressive green filter
- ⌘ Auto and manual focus
- ⌘ Handheld remote
- ⌘ Foot pedal for hands-free image capture



Applications



Hospitals

Forensics



Universities

Telemedicine



Certifications



13458:2012



9001:2008

Mounting Options



Vertical Stand



Swing Arm



Workstation



Laptop



Desktop



Monitors



EHR & EMR Server

Specifications

Camera	SONY high definition camera
Pixels	3.27 mega pixels makes this the highest resolution FDA approved colposcope within its category
Resolution	Compatible with the following screen resolutions: 1920 x 1080, 1280 x 720, 720 x 480, 720 x 576
Image Quality	High definition
Magnification	32X. Optical: 20X, Digital: 12X
Visual Angle	360° camera rotation
Working Distance	7in - 12in
Video Output	HDMI, SDI
Green Filter	Progressive green filter offers a wider range of contrast
Focus	Automatic and manual focus
Image Freeze	Image freeze on screen available with remote, foot pedal, or on camera
Remote	Remote minimizes motion while adjusting the image
Battery	Can work up to 4 hours using battery power



LATERAL
MEDICAL

Contact: (07) 3139 1207 or
sales@lateralmedical.com