

# FIBER NASO-PHARYNGO-LARYNGOSCOPE



*Boundless in Detail!*

# HNO ENT



## FLEXIBLE ENDOSCOPES

### FIBRE NASO-PHARYNGO-LARYNGOSCOPE



Model	FNS-2800	FNS-3400
Insertion Tube	ø 2.8mm	ø 3.4mm
Working Length	300mm	300mm
Fields of View	85°	85°
Angulation Up/Down	180°/180°	180°/180°
Accessories	Leak tester ETO cap, suitcase	Leak tester ETO cap, suitcase
Item No:	FNS-2800G <span style="color: green;">■</span>	FNS-2400G <span style="color: green;">■</span>
	FNS-2800B <span style="color: blue;">■</span>	FNS-3400B <span style="color: blue;">■</span>
	FNS-2900R <span style="color: red;">■</span>	FNS-3400R <span style="color: red;">■</span>

## SPARE PART / ACCESSORIES FOR FIBRE NASO-PHARYNGO-LARYNGOSCOPE

		
Description	Leak tester	PORTLED EMOS mobile light source
Item	200.00040.00	200.02531.00

# HNO ENT

## FLEXIBLE ENDOSCOPE

### SPARE PART / ACCESSORIES FOR FIBRE NASO-PHARYNGO-LARYNGOSCOPE




## LIGHT GUIDE CABLES

### Maxilux pro light guide cable Ø 3,5mm

	1800mm	2300mm	3000mm	Colour
	200.00013.10	200.00014.10	200.00015.10	200.00015.10
200.00013.11	200.00014.11	200.00015.11	200.00015.11	Blue
200.00013.12	200.00014.12	200.00015.12	200.00015.12	Red
200.00013.13	200.00014.13	200.00015.13	200.00015.13	Green
200.00013.14	200.00014.14	200.00015.14	200.00015.14	Black
200.00013.15	200.00014.15	200.00015.15	200.00015.15	Purple
200.00013.16	200.00014.16	200.00015.16	200.00015.16	White
200.00013.17	200.00014.17	200.00015.17	200.00015.17	Orange
200.00013.18	200.00014.18	200.00015.18	200.00015.18	Yellow

### Adapter

	Item No.	Compatible	Description
	200.00003.00	Light Source	Storz® Adapter
	200.02374.00	Light Source	Storz® Adapter
	200.02121.00	Light Source	Fuji® Adapter Xenon
	200.00005.00	Light Source	Wolf® Adapter
	200.02375.00	Light Source	Wolf® Adapter
	200.00006.00	Light Source	Olympus® Adapter
	200.02587.00	Light Source	Olympus® Adapter
	200.00011.00	Light Source	Pentax® Adapter
	200.00007.00	Scope	Storz® / Olympus® Adapter
	200.00008.00	Endoscope	Wolf® Adapter
	200.02588.00	Endoscope	ACMI® Adapter
	200.02704.00	Cable/Endoscope	Storz® / Wolf® Combi



**LATERAL**  
MEDICAL

## CONTACT

**T** +61 7 3139 1207  
**F** +61 7 3036 6187  
**E** [sales@lateralmedical.com](mailto:sales@lateralmedical.com)