

# EPIFAITH®

## Certainty Can Be Effortless



### An epidural space locating syringe with Faith Signal

EpiFaith®, a compact and robust device, assists physicians in locating the epidural space. Intuitively equipped with the Faith Signal visual indicator upon reaching the epidural space, EpiFaith® is beneficial to both junior and senior physicians by transforming the process into a safe and convenient procedure.



### Faith Signal

EpiFaith® provides real-time objective and quantitative pressure information. A visual signal is triggered at the exact moment the needle tip enters the epidural space, indicating a clear endpoint for the advancing of the needle.



### Intuitive Control

By pushing/pulling the plunger to start/stop pressure sensing, EpiFaith® can easily be controlled and monitored by physicians, allowing them to use two hands to advance the needle.



### Feasibility

Equipped with NRFit or Luer connectors, EpiFaith® is compatible with both cutting-edge or mainstream products used in epidural procedures, allowing physicians to maintain their current practices when transferring to EpiFaith®.