


DENTAL LARYNGEAL MASK AIRWAY WITH BITE BLOCK

The **Nasal Laryngeal Mask Airway from Genesis Airway** is the world's first commercially available nasal laryngeal mask. It allows good surgical access for dental, maxillofacial, ear nose and throat, head and neck surgery. A Laryngeal Mask Airway is less invasive than an endotracheal tube (ETT), making it ideal for day surgery.



FEATURES

- Standard flexible Laryngeal Mask Airway.
- Easy clip in/out mechanism allowing repositioning of the Bite Block.
- The Bite Block can be used vertically when waking patient from anaesthesia preventing biting down on Laryngeal Mask Airway tubing.
- Unobstructed surgical access to the entire oral cavity.

DEVICE CODE	DEVICE NAME	SIZE
DFLMA50	Dental Laryngeal Mask Airway	5.0
DFLMA40	Dental Laryngeal Mask Airway	4.0
DFLMA30	Dental Laryngeal Mask Airway	3.0
DFLMA25	Dental Laryngeal Mask Airway	2.5

Units : Minimum Order Quantity 5 and Multiples thereof

- Laryngeal Mask Airway is less invasive than an ETT making it ideal for day surgery. Improved recovery profile allowing an increase in theatre and surgical time efficiency.
- The Laryngeal Mask Airway tubing is locked into the medial side of Bite Block.
- The Bite Block provides a fulcrum over which the Laryngeal Mask Airway may be bent and taped over the cheek, minimizing interference with the surgical team.

