

OBTURATING NASAL ENDOTRACHEAL TUBE

Evidence based studies have shown that the use of an obturating introducer for nasal intubation significantly reduces the risk of incidence and severity of epistaxis to the same, or less than, the incidence when using a red rubber catheter. It also reduced incidence of nasopharyngeal soft tissue damage and nasal pain post intubation.



- The proximal end of the intubation device is manufactured to connect with an endotracheal tube with a posterior facing bevel that allows for easier placement of the endotracheal tube in the trachea without the need for rotation.
- This device makes passage through the nasopharynx a simple technique. If the catheter passes from the nasopharynx to the oropharynx so will the endotracheal tube.

- Soiling of the lumen of the endotracheal tube in the passage through the nasopharynx is prevented as the introducing device obturates the apertures of the endotracheal tube.
- The external diameter of the complex passed through the nose remains the same as the external diameter of the endotracheal tube.
- The Genesis Nasal Obturating ETT is made from softer PVC material, making it easier to use with a bougie. Softer PVC also negates the requirement for thermo-softening prior to nasal intubation.

DEVICE CODE	DEVICE NAME	SIZE
DONETRF80	Obturating Introducer for Nasal Endotracheal Tube - Reinforced Flexible	8.0
DONETRF75	Obturating Introducer for Nasal Endotracheal Tube - Reinforced Flexible	7.5
DONETR70	Obturating Introducer for Nasal Endotracheal Tube - Reinforced Flexible	7.0
DONETRF65	Obturating Introducer for Nasal Endotracheal Tube - Reinforced Flexible	6.5
DONETRF60	Obturating Introducer for Nasal Endotracheal Tube - Reinforced Flexible	6.0
DONETRF55	Obturating Introducer for Nasal Endotracheal Tube - Reinforced Flexible	5.5
DONETRF50	Obturating Introducer for Nasal Endotracheal Tube - Reinforced Flexible	5.0
DONETRAE80	Obturating Introducer for Nasal Endotracheal Tube - Nasal (RAE)	8.0
DONETRAE75	Obturating Introducer for Nasal Endotracheal Tube - Nasal (RAE)	7.5
DONETRAE70	Obturating Introducer for Nasal Endotracheal Tube - Nasal (RAE)	7.0
DONETRAE65	Obturating Introducer for Nasal Endotracheal Tube - Nasal (RAE)	6.5
DONETRAE60	Obturating Introducer for Nasal Endotracheal Tube - Nasal (RAE)	6.0
DONETRAE60	Obturating Introducer for Nasal Endotracheal Tube - Nasal (RAE)	5.5
DONETRAE60	Obturating Introducer for Nasal Endotracheal Tube - Nasal (RAE)	5.0

Units : Minimum Order Quantity 10 and Multiples thereof



**LATERAL
MEDICAL**

1300 144 734
sales@lateralmedical.com